

## Sartorius BCE6202I-1CEU Top-Loading Balance



6200g x 0.01g. From 2021!

Rating: Not Rated Yet

**Price**

Sales price \$985.00

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

### Description

Precision top-loading balance, offering reliable performance for routine weighing tasks with a 6,200 g capacity and 10 mg readability.

### Key Specifications:

Maximum Capacity: 6,200 g

Readability: 10 mg

Repeatability: ? 0.9 mg

Weighing Pan Dimensions: 182 x 182 mm

Adjustment Type: Internal calibration (isoCAL)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)