


## Sartorius MCA36201S-2S00-0 Cubis II High Capacity Balance, 36,200 g x 100 mg



36,200 g x 100 mg

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

### Description

Sartorius Cubis II MCA top-loading precision balances are completely customizable to align with your unique application and compliance requirements while ensuring reliable results. Innovative technology includes automatic motorized leveling, motion controlled and illuminated draft shield, and climate sensor for environmental monitoring to ensure error-free operation. The sizable color touch-screen display allows you to touch, scroll, and swipe with factory-installed weighing applications. The status center shows all information about your balance and environmental conditions as well as detailed help and support in case of warnings and errors.

36,200 g x 100 mg

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

Balances: Sartorius MCA36201S-2S00-0 Cubis II High Capacity Balance, 36,200 g x 100

---