

Sartorius TE 153s



The performance specifications of the Sartorius Talent series of balances set new accents in the class of compact and affordable laboratory balances. Whether you need to weigh in laboratories, vocational training facilities or classrooms, with a Sartorius

Rating: Not Rated Yet

Price

Sales price \$600.00

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

Description

!!!SOLD!!!

The performance specifications of the Sartorius Talent series of balances set new accents in the class of compact and affordable laboratory balances. Whether you need to weigh in laboratories, vocational training facilities or classrooms, with a Sartorius Talent series balance you'll always have just the right equipment.

Capacity	150 g
Readability	.001 g
Reproducibility	+/- .0015g
Linearity	+/- .003 g
Response time	2.5 seconds
PN: with draft shield	TE153S-D
Calibration Weight	100g