

## Savant DNA110 DNA SpeedVac



Are you simply drying down DNA / RNA samples in small volumes of water and ethanol /methanol from 0.4ml to 1.5ml? Then Savant DNA 110 SpeedVac ® System unit may be the perfect choice for you!

Rating: Not Rated Yet

**Price**

Sales price \$1,650.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

### Description

This unit has automatic pump and rotor control. Choice of three drying rates to optimize evaporation. Built in oil-free pump and integrated design miniature design takes up only 1 foot of bench space. Extremely quiet operation. This instrument is capable of using either the 24 place or 36 place microcentrifuge tube rotor. Larger capacity rotor available for 0.5ml tubes. Believe it or not there is even a microplate rotor for the **DNA 110!**

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)