

Savant GP 110-120 Vacuum Pump



The environmentally friendly design of the Savant GP-110 vacuum pump protects users and eliminates hazardous oil waste. Chemically Resistant Design for use with oil pump destroying solvents. Oil-free vacuum pumps create a vacuum without the need for oil,

Rating: Not Rated Yet

Price

Sales price \$750.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

Description

The environmentally friendly design of the Savant GP-110 vacuum pump protects users and eliminates hazardous oil waste. Chemically Resistant Design for use with oil pump destroying solvents. Oil-free vacuum pumps create a vacuum without the need for oil, eliminating the hazards of contaminated oil vapors, oil clean-up and dedicated periodic maintenance of changing oil. Oil-free vacuum pumps are generally much quieter than oil pumps a priority when working within the laboratory environment. They are also much more chemically resistant to vapors passing into the pump. This type of pump is recommended for evaporating samples in acidic environments such as HCl, TFA, Acetic or Formic Acids. Displacement is 36L/min. Maximum vacuum is 7 torr (9mbar)

Instrument Attributes

Additional Specifications

Displacement: 36L / min,
Ultimate vacuum level: 9 mbar,
Dimensions: 33W x 27D x 29H cm,
Weight: 14Kg,
Electrical: 230V / 50-60Hz

Vacuum Pumps: Savant GP 110-120 Vacuum Pump

Biotech Equipment Sales, Inc. based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

[Biotech Equipment Sales, Inc.](#)
226 Miller Ave.
South San Francisco, CA 94088
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com