

Savant Speedvac SPD111V





Deep aluminum TEFLON® coated chamber, accommodates deep-well plates Built in automatic bleeder valve Variable drying rate Connector for vacuum source Time/temperature display incorporating membrane touch pad switches Safety coverlock Manual/au

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

Description

- Deep aluminum TEFLON® coated chamber, accommodates deep-well plates
- Built in automatic bleeder valve
- Variable drying rate
- Connector for vacuum source
- Time/temperature display incorporating membrane touch pad switches
- Safety coverlock
- Manual/auto timed run