

## Savant VLP200 ValuPump



SOLD High Vacuum pumps deliver critical performance for all vacuum requirements. A three position gas ballast control optimizes pump performance to suit all applications.

Rating: Not Rated Yet

**Price**

Sales price \$1,250.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

### Description

**SOLD**

High Vacuum pumps deliver critical performance for all vacuum requirements. A three position gas ballast control optimizes pump performance to suit all applications. The large, sideview sight glass allows for convenient checking of proper oil levels at a glance. The VLP series of pumps is designed to work with all Thermo SpeedVac<sup>®</sup> concentrating and freeze drying systems for maximum solute recovery.

Identical to Edwards RV8