

## Seward Stomacher 3500



Paddle blender with up to 3500ml capacity.

Rating: Not Rated Yet

**Price**

Sales price \$3,450.00

[Ask a question about this product](#)

Manufacturer [Seward](#)

### Description

The Stomacher® 3500 provides a versatile range of tools for large volume blending and extraction where all washing solution or solvent needs to be retained for analysis. This heavy duty machine is built to industrial standards and often used in applications for large volume food microbiology, forensic residue analysis and testing of waterborne parasites.

Sample Volume: 400-3500ml

Dimensions (W x D x H): 19 x 19 x 17.5 in

Power Requirements: 110V, 60Hz, 120W

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**