

## Sheldon SRI3P Refrigerating Incubator



Shel Lab Model SRI3P Refrigerated B.O.D. Incubator 15°C to 30°C (at 20°C ambient) 3 Cu.Ft. (85 L) 120V

Rating: Not Rated Yet

**Price**

Sales price \$1,950.00

[Ask a question about this product](#)

### Description

Shel Lab Model SRI3P Refrigerated B.O.D. Incubator 15°C to 30°C (at 20°C ambient) 3 Cu.Ft. (85 L) 120V

With: Microprocessor Control System; P.I.D. Temperature Control; LED Display of Set point and Chamber Temperature; Interior Electrical Outlet; Fan Assisted/Forced-Air Circulation

### Specifications

Interior Outlet: 110V-1

Agency Approval: TUV-CUE

Shelving: 6 supplied, Adjustable

Temperature Control: 0.1°C Display Increments

Capacity of Standard B.O.D. Bottles: 84

Temperature Uniformity: +/-0.5°C at 20°C

Temperature Range: 15°C to 30°C (at 20°C ambient)

Exterior Dimensions WxDxH in. (cm): 23.5 x 21.0 x 33.7 (59.7 x 53.4 x 85.6)

Interior Dimensions WxDxH in. (cm): 19.0 x 16.2 x 26.2 (48.3 x 41.1 x 66.5)

Capacity: 3 Cu.Ft. (85 L)

## Incubators: Sheldon SRI3P Refrigerating Incubator

---

Nominal Voltage 50/60HZ (V): 100-120/220-240  
Maximum Energy Consumption (W): 360/480  
Door Ajar Alarm: Yes

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)