


## Sigma 4-16KS Refrigerated Benchtop Model w/2-place Deepwell Rotor, 220-240V



Max. capacity [ml]: 4 x 750 ml

Rating: Not Rated Yet

### Price

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [Sigma](#)

### Description

Designed for continuous service, the Sigma 4-16KS is precise and reliable. The centrifuge also features especially good cooling performance and achieves 4°C at maximum speed. The motorized lid lock and the patented user-friendly one-button control (Spincontrol S) make it especially easy to use. A large range of accessories and a variety of rotors made from aluminium for high thermal conductivity allow a broad spectrum of applications.

### Specifications

Max. speed [min<sup>-1</sup>]: 15.000

Max. RCF [x g]: 25.155

Max. capacity [ml]: 4 x 750 ml

Temperature setting range [°C]: -20 – +40

Programs: 60

Acceleration curves: 10 linear, 10 quadratic, 10 user definable

Deceleration curves: 10 linear, 10 quadratic, 10 user definable

Cycle counting: Yes

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)