

## Sorvall F-16/250 Fixed-Angle Rotor



Sorvall F-16/250 Fixed-Angle Rotor for use in Sorvall Evolution RC centrifuge

Rating: Not Rated Yet

**Price**

Sales price \$1,500.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

### Description

#### Rotor F-16/250\*

Capacity (place x ml) - 6 x 250

Max Speed (rpm) - 14,500 in Evolution RC centrifuge

Max.RCF (x g) - 31,940 in Evolution RC centrifuge

Accel/Decel (min:sec) - 2:30/1:45

Tube Size (range, ral) - 3-250

Catalog Number - 11966

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

