

## Sorvall Heraeus 3325B



Sorvall #3325B, Max 13,000 rpm, Max load 24 x 4g, 24 place, 1.5/2mL

Rating: Not Rated Yet

**Price**

Sales price \$150.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

### Description

Sorvall #3325B, Max 13,000 rpm, Max load 24 x 4g

Standard rotor runs 24 tubes in a single row with equalized g-forces for reproducibility, from 1.5 to 2.0 mL tubes to mini-preps and spin columns.

Rotor is for use in Heraeus Pico and Heraeus Fresco centrifuges.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)