

## Sorvall SA-600



Sorvall SA-600 Fixed-Angle Rotor, 12 Places, 50mL each

Rating: Not Rated Yet

**Price**

Sales price \$1,300.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

### Description

SORVALL SA-600 Fixed-Angle Rotor and lid in good condition. For use with Sorvall Evolution RC and RC-5 series centrifuges.

Specifications:

12 places, 50 mL ea.

Diameter: 27cm (10-3/4 in)

Mass (weight): 10.2 kg (22.5 lbs)

Angle: 34 degrees

Number of Places: 12

Maximum Angular Velocity (in rev/min): 16,500/min

Maximum Relative Centrifugal Force (RCF): 39,400

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080

[650-871-5707 office](tel:650-871-5707)

[650-276-7487 fax](tel:650-276-7487)

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)