

Sorvall SM-24 Fixed-Angle Rotor



Sorvall SM-24 Aluminum Fixed-Angle Rotor for use with Sorvall RC-5, RC-5B (inc Plus), RC-6 (inc Plus) and Evolution RC superspeed centrifuges

Rating: Not Rated Yet

Price

Sales price \$850.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

Sorvall SM-24 Aluminum Fixed-Angle Rotor for use with Sorvall RC 6 Plus and Evolution RC superspeed centrifuges.

Specifications -

Capacity 24 x 16mL

Max. RCF 42,717 x g

Angle 28°

Speed 20,500rpm (Outer Row) 51,965 x g (Inner Row)

Diameter (English) Tube 0.71 in.

Diameter (Metric) Tube 18mm

K Factor 434 Outer/591 Inner

For Use With Thermo Scientific™ Sorvall™ RC 6 Plus and Evolution RC superspeed centrifuges

Centrifuge Rotors: Sorvall SM-24 Fixed-Angle Rotor

Type Fixed Angle

Material Aluminum

Length (English) Tube 4.02 in.

Length (Metric) Tube 102mm

Net Weight (English) 17 lb.

Net Weight (Metric) 7.7kg

Radius (English) Max. 4.36Outer3.6 Inner in.

Radius (English) Min. 2.1 in. Outer1.34 in. Inner

Radius (Metric) Max. 110.7mm Outer91mm Inner

Radius (Metric) Min. 53.8mm Outer34.1mm Inner

Min.RCF 25,525xG Outer16,007xG Inner

Biotech Equipment Sales, Inc.

226 Miller Ave.

South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com