

## Stemcell Technologies Easy 250 EasySep Magnet



use with a standard T-75 cm<sup>2</sup> cell culture flask and EasySep™ reagents.

Rating: Not Rated Yet

**Price**

Sales price \$3,450.00

[Ask a question about this product](#)

### Description

Manual cell isolations from large-volume samples, such as leukopaks and whole blood samples, can be scaled up using the Easy 250 EasySep™ Magnet with a standard T-75 cm<sup>2</sup> cell culture flask and EasySep™ reagents.