

Stemcell Technologies Easy 50 EasySep Magnet



Designed to hold a standard 50 mL conical tube, enabling you to process up to 4 x 10⁹ cells in one isolation procedure.

Rating: Not Rated Yet

Price

Sales price \$750.00

[Ask a question about this product](#)

Description

Magnet for column-free immunomagnetic separation

Perform magnetic cell separation procedures easily and efficiently using the Easy 50 EasySep Magnet with EasySep cell isolation reagents. The Easy 50 EasySep Magnet generates a high-gradient magnetic field in the interior cavity that is strong enough to separate cells labeled with magnetic EasySep particles without the use of columns. The Easy 50 EasySep Magnet is designed to hold a standard 50 mL conical tube, enabling you to process up to 4 x 10⁹ cells in one isolation procedure.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

