

Systec V-95



SOLD! Advanced microprocessor controlled toploading autoclave, 16" diameter x 30" deep

Rating: Not Rated Yet

Price

Sales price \$12,000.00

[Ask a question about this product](#)

Manufacturer [Microbiology International](#)

Description

SOLD! Advanced microprocessor controlled toploading autoclave, 16" diameter x 30" deep, accepts up to eight 2-liter flasks, automatic safety door lock, automatic exhaust steam condensation, for liquids & solids, current model, very clean, 3-phase 220V, needs water & air supplies, manufactured in 2004