

## Taitec ECB Electronic Cooling System



Temperature Range +3°C - +20°C (\*) Temp. Control Accuracy Within  $\pm 1^\circ\text{C}$  Standard Attachment Bath Miracle beads bath Temperature Display Digital display Heating-Cooling System Cooling by Peltier device Other Functions Power switch (becomes waiting

Rating: Not Rated Yet

**Price**

Sales price \$525.00

[Ask a question about this product](#)

Manufacturer [Taitec](#)

Description	Temperature Range	+3°C - +20°C (*)
	Temp. Control Accuracy	Within $\pm 1^\circ\text{C}$
	Standard Attachment Bath	Miracle beads bath
	Temperature Display	Digital display
	Heating-Cooling System	Cooling by Peltier device
	Other Functions	Power switch (becomes waiting state)
	Safety device / function	Fuse, overheat protection cutout
	Dimensions	170 x 240 x 220Hmm
	Weight	Body 2.2kg + miracle beads bath 1.8kg
	Power Supply	AC100V 1A
	Exchangeable Sold Bath	Dual side block, hetero block

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotequipmentsales.com](mailto:sales@biotequipmentsales.com)