

Terumo BCT SCD 11B Sterile Tubing Welder - NEW IN PCKG.



Manufactured in 2022. New in packaging.

Rating: Not Rated Yet

Price

Sales price \$13,750.00

[Ask a question about this product](#)

Description

Manufactured: 2022. Model 3NCC986

The Terumo SCD IIB Sterile Tubing Welder, model 3-NCC986, is a closed-system device that can create strong, reliable welds for sterile tubing connections in about 30 seconds. It's compact and portable, and can be easily integrated into a manufacturing process and standard operating procedure (SOP). The SCD IIB is ideal for polyvinyl chloride (PVC) and C-FLEX tubing type. It's also a cleanroom piece of equipment that can weld two C-Flex production hoses in a safe and sterile way.

SCD IIB Welder Connects:

Polyvinyl chloride (PVC) tubing:

Outer diameter: 0.215 inches (5.46 mm)

Wall thickness: 0.033 in 0.84 mm (0.84 mm)

C-Flex tubing:

Outer diameter: 0.25 inches (6.35 mm)

Inner diameter: 0.125 inches (3.18 mm)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com