

## The Viogene miTotal RNA Extraction Miniprep System



The Viogene miTotal RNA Extraction Miniprep System provides a fast method to purify total RNA from various samples such as cells from culture, tissues, whole blood, plasma, serum, biological fluids containing RNA virus etc.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Viogene](#)

### Description

The Viogene miTotal RNA Extraction Miniprep System, which provides a fast method to purify total RNA from various samples such as cells from culture, tissues, whole blood, plasma, serum, biological fluids containing RNA virus etc.. A simple spin-column based method can isolate large RNAs, siRNAs, microRNAs, and viral RNAs without the time-consuming procedure of phenol/chloroform extraction and ethanol precipitation.

### Features and Benefits:

Enhanced small and microRNA yield

Apply all kind of animal cells, tissues, fluids

### Downstream Applications:

Northern blotting

Ploy A+ RNA selection

RT-PCR cDNA synthesis

Transcription profiling

To see Viogene's description of this item: [CLICK HERE!](#)