

## Thermo Arktik Thermal Cycler



96 well x 0.2 ml

Rating: Not Rated Yet

**Price**

Sales price \$1,950.00

[Ask a question about this product](#)

### Description

From individual reactions to high-throughput projects, suitable for a wide range of applications. With easy-to-use graphical interface for protocol setup and preprogrammed protocols for maximum convenience. Manufactured in 2014.

-Block configuration: 96 x 0.2 ml

-Max ramp rate: 3.0°C/s

-Thermal uniformity: +/- 0.4°C at 90°C

-Thermal accuracy: +/- 0.3°C

-Thermal range: 4-99.9°C

-Dimensions: ( W x D x H) 29 x 38 x 29 cm

-Electrical requirements: 100-240V; 50-60Hz

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)