

## Thermo Scientific Model 7450 CryoMed Controlled Rate Freezer



Recent model (2020).

Rating: Not Rated Yet

**Price**

Sales price \$9,500.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

Precise, repeatable freezing results with a temperature ranging from -180 C to 50 C, protecting samples from intracellular freezing. The integrated control system features a freezing chamber, microprocessor control and printer in one benchtop unit, allowing multiple freezing protocols tailored to specific protocols. LN2 tank not included.

### Specifications:

Temp range: -180 - 40C

Capacity: 17L

Access: front

Interior Dimensions: 7 x 12 x 13 in (17.8 x 30.5 x 33 cm)

Exterior Dimensions: 31.7 x 21.2 x 20.4 in (80.5 x 53.8 x 51.8 cm)

Electrical requirements: 12 v / 60Hz / 12.0A 1 Ph

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)