

## Thermo Fisher Scientific -80°C Benchtop Freezer



Ideal for personal or clinical trial applications offering a full 1 cu. ft. (28L) capacity.

Rating: Not Rated Yet

### Price

Sales price \$3,750.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

### Description

ULT Freezer, 1CF(28L), 115V/60Hz; ideal for personal or clinical trial applications offering a full 1 cu. ft. (28L) capacity.

- Dimensions (L x W x H) Exterior: 25.5 x 30.75 x 20 in. (64.8 x 78.1 x 50.8cm)
- Dimensions (D x W x H) Interior: 12 x 15.3 x 9.6 in. (30.5 x 38.7 x 24.5cm)
- Plug Type: NEMA 5-20
- Temperature Range: ?50° to ?80°C
- Footprint (English): 0.106 cu. ft.
- Amperage: 16/20
- Display: Push Button LED
- Shipping Weight (English): 190 lb.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)