

Thermo Fisher Scientific TSX 2320FA High-Performance -20°C Manual Defrost Freezer





Dimensions (D x W x H): 28.5 x 24 x 58 in.

Rating: Not Rated Yet

Price

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Plug type: NEMA 5-15

Electrical requirements: 115V, 60Hz

Dimensions (D x W x H): 28.5 x 24 x 58 in

Unique V-drive technology is designed to detect usage patterns such as door openings when a higher compressor speed is needed and periods of stability where the compressor runs at a lower speed, saving energy without compromising protection. Cold wall convection cooling for uniform and stable temperatures. Easy to use microprocessor control with interior light, setpoint security and alarms. Standard door locks provide comfort and assurance. Ergonomically designed door handle with a soft-touch grip. The controller is angled to help prevent neck strain. Automatic pressure equalization port that enables fast reentry. Easy to clean and adjust epoxy coated shelves. Whisper quiet with a sound level of 52 dBA, enabling end users to work alongside the refrigerator without distraction. Self-closing door with 90° stop to assist with inventory loads.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com