

## Thermo Forma VLP 80 Vacuum Pump



SOLD Three-position gas ballast control optimizes pump performance to suit all applications Large, sideview sight glass allows for convenient checking of proper oil levels at a glance All Rotary Vane Oil Pumps come standard with: Muffler (OM

Rating: Not Rated Yet

**Price**

Sales price \$850.00

[Ask a question about this product](#)

Manufacturer [Thermo Forma](#)

### Description

**SOLD**

- Three-position gas ballast control optimizes pump performance to suit all applications
- Large, sideview sight glass allows for convenient checking of proper oil levels at a glance
- **All Rotary Vane Oil Pumps come standard with:**
  - Muffler (OME190 Oil Mist Eliminator)
  - Two Oil Filters (MF190)
  - Quick Drain Kit