

Thermo IEC Centra CL2 General-purpose benchtop centrifuge with 4-place swinging bucket rotor



Thermo IEC Centra CL2 General-purpose benchtop centrifuge with 4-place swinging bucket rotor The IEC Centra CL2 general-purpose benchtop centrifuge incorporates an impressive array of features that ensure safe, convenient operation and consistent results for your low-speed applications.

Rating: Not Rated Yet

Price

Sales price \$1,100.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Thermo IEC Centra CL2 General-purpose benchtop centrifuge with 4-place swinging bucket rotor

The IEC Centra CL2 general-purpose benchtop centrifuge incorporates an impressive array of features that ensure safe, convenient operation and consistent results for your low-speed applications.

- 4-place swinging bucket rotor
- Impressive range of speeds – up to 4,770 x g
- Small footprint requires minimal bench space
- Lightweight and easy to transport
- Sealed Aerocarrier™ adapters ensure biocontainment
- High capacity – Spin up to 28 x 5/7 mL (13 mm) blood tubes and 8 x 50 mL conicals
- User-friendly controls allow rapid time and speed adjustments
- Microprocessor-based controls ensure consistent, accurate results