

Thermo Labsystems Multidrop Micro



1-50 μ l into 384- and 96-well microplates; stand-alone instrument that can be integrated easily with robotic systems.

Rating: Not Rated Yet

Price

Sales price \$3,500.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Using specifically developed technology Multidrop Micro rapidly and accurately dispenses micro-volumes of 1-50 μ l into 384- and 96-well microplates. Using increments of just 1 μ l Multidrop Micro provides a new level of precision for your dispensing needs. It is easy to select the exact volumes you require with the convenient onboard switches. This stand-alone instrument has the flexibility to be integrated easily with robotic systems.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com