


## Thermo M-10 Titerplate Rotor for ST8 Centrifuges



Capacity: 4 standard microplates or 2 deep-well (midi) microplates; Maximum Speed: 4,400 rpm

Rating: Not Rated Yet

**Price**

Sales price \$1,150.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

The M-10 microplate rotor is a swinging-bucket rotor engineered for use in Sorvall ST8, ST8R, Megafuge 8/8R, and SL8/SL8R centrifuges. It supports high-throughput microplate workflows, handling up to 4 standard microplates or 2 deep-well plates with interchangeable plate carriers. The rotor uses Thermo's Auto-Lock system for quick, tool-free installation, and optionally supports biocontainment with ClickSeal lids.

### Specifications

Capacity: 4 standard microplates or 2 deep-well (midi) microplates

Maximum Speed: 4,400 rpm

Maximum Relative Centrifugal Force (RCF): 2,576 × g

Rotor Radius: 119 mm

K-Factor: 13,103

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)