

Thermo Masterflex Peristaltic Pump 1300-3100



With Masterflex L/S EasyLoad Head II model 77201-60

Rating: Not Rated Yet

Price

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

0.02-100 rpm 0.1 HP; 115/230 AC, 50/60 hz, 1 Ph, 2.2/1.1 A

With the 77201-60 Masterflex L/S EasyLoad Head II:

Pump series: L/S

Number of channels: 1

Max flow rate: 2900 mL/min

Min flow rate: 0.034 mL/min

Max flow rate per channel: 2900 mL/min

Min flow rate per channel: 0.034 mL/min

Volume per revolution with L/S 13 tubing: 0.06 mL

Volume per revolution with L/S 14 tubing: 0.21 mL

Volume per revolution with L/S 16 tubing: 0.8 mL

Volume per revolution with L/S 25 tubing: 1.7 mL

Volume per revolution with L/S 17 tubing: 2.8 mL

Volume per revolution with L/S 18 tubing: 3.8 mL

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com