

## Thermo NapCOIL UF 500 -80 Freezer



The NapCOIL\*, cooling-coils-in-the-shelf design is a new concept in this field. With the cooling coils in direct contact with the sample, freezing is very rapid, and sample temperature uniformity is improved.

Rating: Not Rated Yet

### Price

Sales price \$4,850.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

The NapCOIL\*, cooling-coils-in-the-shelf design is a new concept in this field. With the cooling coils in direct contact with the sample, freezing is very rapid, and sample temperature uniformity is improved.

- . Preformed insulation blocks eliminate the possibility of uneven insulation and "hot spots"
- . Inner doors seal on the shelf edge so that cold air will stay in unopened compartments
- . Two oversized 3/4 hp compressors and an oversized first stage condenser eliminate the need for auxiliary cooling fans
- . Robust design assures long trouble-free service, lower power requirements and less noise
- . Temperature range: -65° to -85°C (-85° to -121°F)
- . Four shelves