

## Thermo Savant Discovery SpeedVac



Thermo Savant Discovery large volume SpeedVac w/multi-tiered rotor

Rating: Not Rated Yet

**Price**

Sales price \$18,500.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

Description

**!!! SOLD !!!**

This Thermo-Savant Discovery SpeedVac is a high capacity, high throughput centrifugal concentrator designed for combinatorial chemistry and its associated solvents. High quality components handle aggressive solvents such as TFA, Acetonitrile, DMSO, Methanol and Methylene Chloride. Interchangeable multi-level rotor can hold up to 120 shallow well microtiter plates or whole test tube racks in a single run. The compact design allows the complete unit to fit inside a conventional walk-in fume hood without modifications. A built-in 100 L/min oil-free, maintenance-free vacuum pump and a 7.5L capacity cold trap make it a complete, stand-alone system with plug-in and use convenience.

**Please call us at (650) 871-5707 to check availability of this product!**