

Thermo Savant SPD2030 Integrated SpeedVac Concentrator System



Manufactured in 2023.

Rating: Not Rated Yet

Price

Sales price \$14,950.00

[Ask a question about this product](#)

Manufacturer [Thermo Savant](#)

Description

The SPD2030 SpeedVac concentrator is a complete system that integrates an oil-free vacuum pump, a refrigerated vapor trap and high-capacity centrifugal concentrator into a single compact package. Ideal for solvent evaporation, sample concentration and drying. Features an advanced front panel with dual timer for automatic or manual control of run conditions along with preset and user-defined program storage.

Specifications

Chamber Temperature Range: 35°C to 80°C

Refrigerated Trap: -50°C

Maximum Vacuum: <10 Torr

Included Components: Glass vapor trap and RH50-28-60 rotor (supports up to 50 x 28x60 mm scintillation vials)

Electrical Requirements: 220V, 60Hz, 8A

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080

650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com