

Thermo Scientific Heratherm Oven, 2.2 cu ft.





Temperature Range (Metric): 50° to 330°C

Rating: Not Rated Yet

Price

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Model OMH-60S

2 shelves

Effectively protects samples with an extra layer of security for applications where process reliability and security are paramount. A standard over-temperature alarm and an additional under-temperature alarm ensure that your sample stability. A door lock prevents disruption or tampering of samples, while a door alarm notifies the operator when the door is left open accidentally.

Features

Broad temperatures range—up to 330°C.

Advanced digital timer for daily or weekly on / off cycles.

Stainless-steel interior is easy to clean and corrosion resistant (steel quality AISI 304).

Programmable controller for temperature ramps and dwells.

Up to 10 programs can be saved, each with 10 discrete steps.

Electronically controlled fan speed and damper position can be incorporated.

Programs can be repeated automatically.

Access port allows the introduction of sensors for independent data monitoring.

Simple calibration routine ensures temperature accuracy over time.

Boost function enables rapid heating —no need to operate the oven 24/7.

Specifications

Cabinet Material: Coated steel

Capacity (Metric): 62L

Dimensions (D x W x H) Interior: 13.5 x 12.9 x 20 in. (34.3 x 32.8 x 50.8cm)

Dimensions (L x W x H) Exterior: 22.2 x 20.9 x 28.3 in. (56.5 x 53 x 72cm)

Electrical: Amperage: 12A; Hertz: 60Hz; Voltage: 120V; Wattage: 1440w

Shelves: 2 supplied/13 max

Temperature Range (Metric): 50° to 330°C

Temperature Control: Programmable Microprocessor control with vacuum fluorescent display

Timer: weekly / real time / hour

Temperature Stability: ±0.2°C

Type: Mechanical convection

Plug Type Nema 5-15

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax
sales@biotechequipmentsales.com