

Thermo Scientific MaxQ 2000 CO2 Digital Shaker



Small footprint. Variable speed control from 15 - 500 rpm

Rating: Not Rated Yet

Price

Sales price \$1,650.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

CO2-resistant shaker designed for environments containing up to 10% carbon dioxide.

- Variable speed control from 15 - 500 rpm
- Speed and time are controlled on individual LED displays
- Shaking speed ramps up slowly to prevent accidents and sudden spills
- Brushless DC motor does not require maintenance and provides long-life under heavy use
- Shaker discontinues operation and emits an audible alarm if platform is loaded unevenly or the speed fluctuates erratically
- Continuous or timed operation (in minutes or hours)
- Load capacity up to 35 lbs. (15.9 kg)
- Small footprint conserves bench space
- 19 mm (0.75") orbital drive
- Can be recalibrated using a digital tachometer
- Includes RS232 port for PC connection
- Requires a dedicated platform (with predrilled clamps) or a universal platform with individual clamps
- Options include larger platforms, stacked platforms, and clamps for various vessels

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.

South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com