

Thermo Scientific Precision GP 05 Water Bath



5L; temp - ambient to 100C

Rating: Not Rated Yet

Price

Sales price \$400.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

General purpose water bath -- designed to maintain water temperature from ambient to 100°C. Over-temperature safety circuitry is designed to prevent thermal runaway, while new auto-on and auto-off timers allow you to optimize operation schedules. Benefit from outstanding chemical and corrosion resistance with epoxy powder-coated exterior, and easily clean the chamber with its seamless stainless steel interior.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com