

## Thermo Scientific Savant SpeedVac SPD1030 Integrated Vacuum Concentrator



With a 1.5 mL microcentrifuge-tube rotor.

Rating: Not Rated Yet

### Price

Sales price \$7,750.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

The Thermo Scientific Savant SpeedVac SPD1030 is an integrated vacuum concentrator designed for rapid drying and concentration of aqueous and compatible organic-solvent samples, including DNA/RNA, proteins, enzymes, and HPLC fractions. This unit is equipped with a 1.5 mL microcentrifuge-tube rotor.

### Key Features/Specifications

Rotor: 1.5–2.0 mL microcentrifuge-tube rotor; the RH40-11 version holds 40 tubes with approximately 1.2–1.6 mL working volume per tube.

Temperature range: Ambient to 80°C.

Maximum vacuum: <10 Torr (13 mbar).

Integrated system: Includes a diaphragm vacuum pump, –50°C cold trap, and vacuum gauge, eliminating the need for separate external components.

Versatile capacity: The SPD1030 can accommodate up to 120 × 1.5 mL tubes, depending on the rotor installed, as well as 15 mL and 50 mL tubes and microwell plates with appropriate rotors.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

