

## Thermo Scientific SLA 1500 Rotor Super-Lite



Fixed-angle rotor; 15,000 RPM.

Rating: Not Rated Yet

**Price**

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

### Description

Aluminum 6 x 250 ml fixed-angle rotor with a maximum speed of 15,000 rpm; autoclavable at 121°C. For use with Sorvall RC-6 Plus series, Thermo Scientific Sorvall Evolution RC series, and Sorvall RC-5 Plus series centrifuges (including the RC-5, RC-5B, RC-5B Plus, RC-5C, RC-5C Plus, and RC-26 Plus).

[Biotech Equipment Sales, Inc.](#)

226 Miller Ave.

South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)