

Thermo Scientific Standard Tube Dispensing Cassette, New in Box!





For Multidrop Combi, Multidrop 384, Multidrop DW.

Rating: Not Rated Yet

Price

Sales price \$425.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

This standard tube dispensing cassette is designed for use with the Multidrop Combi+ and Multidrop Combi SMART+ dispensers. It is suitable for dispensing volumes ranging from 5 μ L up to 2,500 μ L, with a focus on volumes of at least 20 μ L. The wide tubing allows for high-speed dispensing of large volumes, making it ideal for applications requiring precise and efficient liquid handling.

Specifications

Product Type: Tube Dispensing Cassette

Autoclavable: Autoclavable (50 times at 121°C, 1 bar, 20 min.)

Description: Standard tube dispensing cassette

Material: Polypropylene Tip, Silicone Tubing

Volume (Metric): 5 to 2500 μ L

For Use With (Equipment): Multidrop Combi, Multidrop 384, Multidrop DW

Volume (Metric) Dead: <7 mL

Diameter (Metric) Tip: 0.5 mm; Tubing: 1.3 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com