

Thermo Scientific Thermal Mixer



w/ 24 x 1.5 mL block

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Fast heat-up and cooling times and quick to reach specified mixing speed

Temperature control range: 15°C below ambient to 100°C

Temperature setting range: +4°C to 100°C

Shaking speed: 250 to 1400 rpm

Temperature stability: $\pm 0.1^\circ\text{C}$

Even mixing across the thermoshaker block

Excellent uniformity across the block

Digital timer with sound alarm: 1 min to 96 hours

Temperature calibration function - compensate for differences in the thermal behavior of different tubes and samples

Thermal Cyclers: Thermo Scientific Thermal Mixer

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax

sales@biotechequipmentsales.com