

## Thermo Scientific TSHP series FMS High-Performance Lab Refrigerator



23 cu. ft.; 651, Flammable Materials Storage

Rating: Not Rated Yet

**Price**

Sales price \$3,850.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

### Description

Model number: TSFMS2305A

23 cubic feet, 115 volts, 60 hertz, 13.5A

Protect flammable materials, solvents and intermediates using Thermo Scientific™ TSHP FMS High-Performance Laboratory Refrigerators, which meet stringent U.S. Pharmacopeia standards for drug storage as well as the NFPA safety provisions for the storage of flammable liquid

### Specifications:

Capacity: 23 cu. ft.; 651L

Temperature Range: 1° to 8°C

Shelves: 3

Shipping Weight (English): 425 lb.

Defrost: Manual

Door Style: Solid

Dimensions (Ext.) (H x W x L): 78.6 in. x 28 in x 37 in

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)