

Thermo TSX40086A Ultralow -80C Freezer, 115V, 19.4 cu/ft, from 2018



TSX Series Ultra-Low Freezer, 400-Box Capacity, 115V, 60Hz

Rating: Not Rated Yet

Price

Sales price \$6,750.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

TSX Series Ultra-Low Freezer, 400-Box Capacity, 19.4 cu ft.; 115V, 60Hz

Specifications:

Capacity: 19.4 cu. ft. / 549L; 165 lb / 75kg (shelf)

Dimensions Exterior (L x W x H): 37.6 x 28.3 x 78 in. (95.5 x 71.9 x 198.1cm)
Dimensions (D x W x H) Interior: 28.3 x 23.1 x 51.2 in. (71.9 x 58.8 x 130.1cm)
Plug Type: NEMA 5-15
Temperature Range: -50° to -86°C
Type: Ultra Low Freezer, Upright
Electrical Requirements: 115V, 60Hz
Smart View Compatibility
Weight: 573 lb. / 260kg
Holds: 400 2-inch boxes
2mL Vial Capacity: 40,000
Insulation: Water-Blown Foam (SNAP-compliant)
Refrigerant: R290 in First Stage, R170 With R290 Mix in Second Stage (SNAP compliant)
Shelves: 3
Shipping Weight: 665 lb. / 302kg
Defrost: Manual
Door Style: 1 outer, 2 inner

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com