

ThermoSavant SPD131DDA System



This ThermoSavant SpeedVac system, consisting of the SPD131DDA concentrator, the RTV 4104 refrigerated trap, and the OFP400 teflon oil-free vacuum pump, is solvent resistant and ideal for drug discovery applications. All 3 components were manufactured in

Rating: Not Rated Yet

Price

Sales price \$9,750.00

[Ask a question about this product](#)

Manufacturer [Thermo Savant](#)

Description

!!! SOLD !!!

This ThermoSavant SpeedVac system is solvent resistant and ideal for drug discovery applications. All 3 components were manufactured in 2006.

The SPD131DDA SpeedVac concentrator is microprocessor controlled, with digital displays for time, temp, and vacuum. Has a teflon coated interior w/stainless steel upper magnetic assembly, radiant cover & heated lid. It's Resistant to TFA, DMSO, and other solvents used in combinatorial chemistry applications.

The RVT4104 trap gets down to -104C for trapping low-boiling solvents such as Chloroform or Methylene Chloride and has a digital temp display.

The OFP400 teflon oil-free vacuum pump is recommended for evaporating samples in acidic environments such as HCl, TFA, Acetic or Formic Acid, aggressive/volatile organics (methylene chloride, acetone, hexane) and strong bases (ammonium hydroxide).

This system is located at our facility in South San Francisco & comes w/our 90-day warranty.

Biotech Equipment Sales, Inc. is a new company, based in the Biotech capital of the world, dedicated to bringing high quality pre-owned laboratory instrumentation to our scientific community.

Please call us at (650) 871-5707 to check availability of this product!