

ThermoScientific 4 Point Magnetic Stirrer



Thermo Scientific Stirrer Magnetic 4 Points Biosystem 4 Direct **New in Box**

Rating: Not Rated Yet

Price

Sales price \$1,000.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Thermo Scientific VarioMag Biosystem 4 - stirrers generate little to no heat. Thermally-sensitive cultures and cell suspensions can gently be stirred with complete confidence.

Biosystem 4 Direct stepper-drive units turn stirring bars gently and evenly to preserve valuable cultures. The Biosystem 4 Direct offers 4 synchronized stirring points in a rectangular arrangement.

- Sealed stainless steel housing for use in high humidity environments
- Compact, flat design
- Automatic start-up for safe stirrer acceleration
- Easy to clean under running water

Specifications -

Model: Biosystem 4 Direct

Stirring Points: 4

Gap Between Points (mm): 185

Stirring Volume per Point (mL): 25-2,000

Stirring Power Setting (W): 4 x 0.5

Shakers/Stirrers/Hotplates: ThermoScientific 4 Point Magnetic Stirrer

Speed (rpm): 5-120

Operating Conditions(C): -10 to [+40@95%rH](#)

Protection Class: IP 64

Dimensions W x D x H (mm): 330 x 330 x 60

Weight (kg): 11.0

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com