

Tri-Philizer MP Tri15C0T001/MPPR/VP400C



Pilot Plant FT TriPhilizer has three tray chambers with 5 temperature-controlled shelves in each chamber with capacity to freeze dry 15L simultaneously. The condenser temperature is -85°C .

Rating: Not Rated Yet

Price

Sales price \$12,500.00

[Ask a question about this product](#)

Manufacturer [Pilot Plant](#)

Description

!!!! - RECENTLY SOLD - !!!!

Pilot Plant FT TriPhilizer has three tray chambers with 5 temperature-controlled shelves in each chamber with capacity to freeze dry 15L simultaneously. The condenser temperature is -85°C . Applications:

- solvent removal;
- preservation of bacteria or other biologicals;
- stabilization of bulk chemicals;
- protein preservation.