

## Turner Designs TD-700 Fluorometer



The TD-700 Laboratory Fluorometer is a compact, low-cost benchtop fluorometer that features a unique multi-assay filter cylinder, wide wavelength range, menu-driven software, and multi-point calibration.

Rating: Not Rated Yet

**Price**

Sales price \$2,000.00

[Ask a question about this product](#)

Manufacturer [Turner Design](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

- Multi-Assay Filter Cylinder houses eight easily-changeable filters for four different applications.
- Wide wavelength range: measures fluorescence from UV to red range
- Automatic Range Optimization: range, gain, and resolution are set automatically - eliminates operator error from working in the wrong range.
- Multi-point calibration: stores a 1-to-5 point calibration curve for quantitative determinations.
- Kinetics Software: outputs data to a computer or printer at selectable intervals.
- Versatile Sample Compartment accommodates 25 x 150 mm round, 13 x 100 mm round, 10 x 10 mm square cuvettes, 100 µl minicell, or 9 µl microcell.
- Menu-Driven Software guides the user through calibration and analysis.
- RS-232 for serial data output to a printer or computer.