

## Vacuubrand BVC Professional Fluid Aspiration System



Multiples in stock.

Rating: Not Rated Yet

**Price**

Sales price \$1,375.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

### Description

**Features:**

Powerful and long-lasting chemistry diaphragm pump for universal economical use

Extremely compact, space saving and low noise and thus the ideal tool in the workplace

Easy operation due to electronically adjustable suction power via touch panel, for strong or sensitive suction

Autoclavable collection bottle designed to reduce aerosol and foam formation; with hydrophobic protection filter 0.2 µm, autoclavable quick-couplings available as accessory

Ready to accept a second aspiration hand controller (VHCpro) to support two users

**Specifications:**

## Vacuum Pumps: Vacuubrand BVC Professional Fluid Aspiration System

---

Max. pumping speed 50 Hz: 0.7 m<sup>3</sup>/h (Vacuum pump)  
Max. pumping speed 60 Hz: 0.5 cfm (Vacuum pump)  
Ultimate vacuum: 150 mbar / 112 torr  
Vacuum controller: Touchpanel  
Control range: 150 - 850 mbar (adjustable)  
Collection bottle: 4l Polypropylene  
Max. back pressure (abs.): 1.1 bar  
Inlet connection: Suction controller  
Outlet connection: Hose nozzle DN 8 mm / silencer  
Rated motor power: 0.04 kW  
Rated motor speed at 50/60 Hz: 1500/1800 min<sup>-1</sup>  
Protection class: IEC 60529IP 40  
Weight: 7.7 kg  
Noise (sound pressure level) in dBA at 50 Hz or 1500 rpm/62% (VARIO)/1500 rpm (VARIO-SP)/12500 Upm (VACUU·PURE®): 49 dBA  
NRTL certification: No  
Rated mains voltage range 1: 15 V  
Mains frequency 1: 50-60 Hz

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)