


Vacuubrand Chemistry Hybrid Pump RC 6



Made in 2016. Max. pumping speed: 5.9/6.9 m³/h

Rating: Not Rated Yet

Price

Sales price \$2,500.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

Description

Max. pumping speed: 5.9/6.9 m³/h

Reduced internal corrosion, even when working with corrosive vapors
Oil changes typically reduced 90% or more compared with rotary vane pumps alone
Most economical solution: In practical operation a cold trap is often no longer necessary
Ease of maintenance due to telescopic design

Ultimate vacuum: 2×10^{-3} mbar
Ultimate vacuum with gas ballast: 1×10^{-2} mbar
Max. pumping speed 50 Hz: 5.9 m³/h

See serial plate for more information.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com