

Vacuubrand MD 4C Vacuum Pump System



SOLD Applications include evacuation, evaporation, and pumping of gasses and vapors in chemical, biological, and pharmaceutical laboratories.

Rating: Not Rated Yet

Price

Sales price \$3,950.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

Description

SOLD

Ideal for high vacuum requirements with high boiling solvents; typical applications are rotary evaporators and drying ovens. The separator at the inlet (AK) retains particles and liquid droplets; the waste vapor condenser at the outlet (EK) is highly efficient and compact. The condenser enables near 100 percent solvent recovery, 100-115V, 50/60Hz, 120V, 5.7A.