

Vacuubrand MZ 2C Vacuum Pump



Vacuum pump with Integrated Controller and Condenser.

Rating: Not Rated Yet

Price

Sales price \$2,950.00

[Ask a question about this product](#)

Manufacturer [Vacuubrand](#)

Description

Vacuubrand MZ 2C Vacuum Pump w/ Integrated Controller and Condensor -

Vacuubrand Vario oil-free vacuum systems automatically find and follow boiling points for minimum oversight, greatly reduced bumping, and shorter process times in rotary evaporation, centrifugal concentration, and other evaporative applications.

All systems feature a controller, diaphragm vacuum pump, and solvent recovery accessories (inlet catchpot and outlet condenser). Controller continually adjusts pump speed to changing vapor conditions, for optimum vacuum with no preprogramming or hysteresis. Vacuum depth is displayed in operator's choice of torr, mbar, or hPa.

Integral diaphragm pumps utilize planar diaphragms to provide high flow rates in working vacuum range. Pumps are equipped with PTFE flowpaths for excellent chemical compatibility and long unit life.

Biotech Equipment Sales, Inc. based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

[Biotech Equipment Sales, Inc.](#)
226 Miller Ave.
South San Francisco, CA 94080
650-971-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com