

## Virtis FreezeMobile 25EL



Virtis FreezeMobile 25EL large volume freeze dryer on wheels, to -85C SOLD!!!

Rating: Not Rated Yet

**Price**

Sales price \$12,750.00

[Ask a question about this product](#)

Manufacturer [Virtis](#)

Description

**!!! SOLD !!!**

This Virtis FreezeMobile 25EL large volume freeze dryer on wheels, to -85C comes w/a filter trap, uses one button start up of systems in proper sequence & hot gas defrost that allows easy ice removal & rapid turnaround time. Power save function reduces power consumption & extends vacuum pump life. 208/230V, 60 Hz includes Trivac D16B vacuum pump & two-tiered manifold w/spigots.

SOLD!!!

**Please call us at (650) 871-5707 to check availability of this product!**

